TREKKING MBH
Total Knee Arthroplasty

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Total prosthetic knee system

The Trekking total prosthetic knee system is an integrated system, providing the surgeon a wide range of options to choose from the appropriate solution for each patient and pathology, thus meeting the requirements of each single clinical case.

Trekking is the outcome of a 7-year European multi-center research programme that included hundreds of clinical tests. Actually, Trekking is an high-level system, only comparable with similar products manufactured by the best biomedical companies in the world.

The evolution in the studies regarding knee replacement systems and the updated and proven design of components lead to an extended implant longevity, reproducing a construction of the prosthetic joint as similar as possible to the natural one, with very good performances under the biomechanic, cinematic and articular aspects.

The Trekking System, thanks to its features of biomechanic implant, avails itself of the different fixation opportunities, i.e. the mechanical primary anchorage and the biological secondary.

By Trekking Systems you can obtain a perfect identification of these parameters; an effective means to correct limb dysmetria and as well as to allow The Trekking System, thanks to its precise and easy instrument set, allows an exact identification of the anatomic parameters, granting an effective correction of dysmetria and a better muscular functionality, helping surgeons to conform to the parameters he fixed during the pre-operative planning and maximizing the implant's performance.

Sophisticated manufacturing processes and attentive and particularized control procedures guarantee for the total and full reliability of the products.

Trekking knee system includes:

- three different femoral components - CR, PS and SPS Revision, with the possibility of adding stabilizing stem extension, offset adapters, as well as augments;
- four different tibial components, with inner cone for UHMWPE rotating mobile bearings' peg, rotating and rotating-translating bearings, fixed bearings and Revision fixed bearings, and options for offset adapters, as well as augments;
- related tibial inserts, and further constrained reinforced inserts for coupling with SPS Revision femoral component and fixed tibial plate;
- patellar components;
- two sets of stabilizing, interchangeable titanium stem extensions for Revision femoral and tibial component: one for cemented implantation, mirror polished and with a special design for avoiding problems in cementing, and one with roughened surface and an appropriate design for achieving a correct grip where cementless fixation is indicated;
- a set of augments in different sizes and forms (posterior and distal femoral, block and wedge tibial augmentations);
- titanium offset adapters;
- CR and PS femoral components and tibial plates for fixed bearings with pure titanium TiCoat coating obtained with plasma spray technique for a better osteointegration, suited for cementless implants.

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Total Knee Replacement System

The TREKKING MBH System features a rotating insert in a mirror-polished tibial plate and an articular surface perfectly congruent with the corresponding surface on the femoral component. The mobile bearing knee system has three main advantages:

- Reduced polyethylene wear, thanks to an optimised articular contact area.
- Improved implant function: each movement is accommodated by a different joint. Moreover, each surface has been specifically designed for a dedicated movement, thus considerably improving both wear performance and joint function.
- Tibial plate position is a less critical parameter, since optimal alignment can be naturally reached by the bearing.

The mobile bearing knee system is indicated in relatively young and active patients with good ligaments.

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The TREKKING CR cemented femoral component is manufactured in a CoCrMo alloy. It is an ACL-retaining system and therefore indicated in patients with ligaments in good conditions.

Components are available in cemented and uncemented versions with a VPS treated (Vacuum Plasma Spray) TiCoat.

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
The TREKKING PS cemented femoral component is manufactured in a CoCrMo alloy. It is an ACL-sacrificing system and is provided, with a central cam to ensure natural femoral roll-back during gait starting from 40° of flexion.

Components are available in cemented and uncemented versions with a VPS treated (Vacuum Plasma Spray) TiCoat.

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Circular sagittal profile

The contact area of the Trekking condyle has been designed with a single radius of curvature in the most frequently loaded zone (70° in flexion and 10° hyperextension).

Any change in the radius of curvature in this zone leads to a reduction in the contact area between the condyle and the PE insert, thus increasing the risk of insert wear. It would also produce an unbalanced ligament tension in the transit area between flexion and extension, thus making simple daily actions such as standing up from a chair or pressing the car pedals less comfortable.

In the TREKKING system, the circular profile rationale is extended to include also the frontal section.

In fact, frequent oscillations occur (lift-off) also along this plane, during which articular surface contact must be perfectly maintained, in order to prevent early polyethylene wear.

Therefore, condyles have been designed with a radius of curvature identical to the sagittal profile, so that each condyle is a section of a sphere.

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
Bearing

The TREKKING CR insert is made of Ultra high Molecular Weight Polyethylene (UHMWPE), sterilized with ethylene oxide. Similarly to the corresponding femoral component, it is anterior cruciate-retaining and therefore indicated in patients with good ligaments.

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
Bearing

The TREKKING PS insert is made of Ultra high Molecular Weight Polyethylene (UHMWPE). Similarly to the corresponding femoral component, it is an anterior cruciate-sacrificing system provided with a central cam for natural femoral roll-back during gait starting from 40° of flexion.

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The TREKKING CR insert has been designed with chamfers, both anteriorly to prevent patellar tendon impingement, and posteriorly to improve Range Of Motion (ROM).


Thanks to these features and a particularly accurate femoral component design, the TREKKING system makes it possible to achieve a mechanical ROM (i.e. with no interference of soft tissues) of up to 150°.

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
Tibial component

The TREKKING cemented tibial component is manufactured in a CoCrMo alloy by micro fusion and machining. The plate is mirror-polished to reduce back-side wear and fitted with a central hole to accommodate the UHMWPE insert peg. For better anatomical congruency, the keel comes in five different sizes. Posterior plate tilt is 0° and allows for better insert mobility.

The TREKKING MBH tibial plate system includes a 3 mm plate for considerable bone sparing. Moreover, a finite element method (FEM) has been used to design the keel in such a way as to ensure optimal mechanical strength.

Components are available in cemented and uncemented versions with a VPS treated (Vacuum Plasma Spray) TiCoat.

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
Surgical sequence

This surgical sequence has been studied to obtain at once both the maximum precision and an easy recovering of any unsatisfactory situations.

The tibial resection gives the surgeon a good benchmark to perform the femoral cuts, and more room to operate in the femoral part. To leave the chamfers as the last cuts, allows a much more precise measurement of the articular gap, and a much easier recutting, if any is needed.

A functional system of augmentations of the trial balance helps the surgeon to select the correct thickness of the insert.

At any stage, the surgeon may verify the articular alignment with a metal rod inserted in the proper eyelets present in most of the instrumentation parts.

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

Main Advantages

- It is possible to check the flexion and extension balancing before chamfering the femur.
- Any correction done at this stage does not compromise the bone-implant fit.
- More room to operate in the femoral part
- Good benchmark to perform the femoral cuts
- Reduced number of passages and devices (steps 1b and 3a are made with single masks)

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The proper height of the tibial cut is set by mean of the two-tips Tibial Stilus, which has to be placed on the lower healthiest part of the tibial plate. One tip is marked as "Slot" and must be used if the surgeon intends to cut using the fissure of the cutting block; the other one is marked as "Top", and must be used if a cut in the superior surface is going to be performed.

When the Tibial Cutting Block is in the correct position, two nails must be placed in the "zero" holes. The "+2" and "+4" holes are for recutting. The universal Drill Adapter (1) may be used. The setup of the tibial structure may be verified by means of the Resection Tester. The cutting block can be freed by the simple turn of the frontal lever (1).

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Comfortable push-buttons allow the orientation in the proximal-distal, medio-lateral and sagittal plans.

For a better stability, after the alignment system is removed, a nail may be inserted in the hole in the upper part of the fixation arm.

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This is the SAMO TREKKING femoral system, one of the most advanced concept for femoral measurement, alignment and cut. The system allows the surgeon to perform all the measurement and the cuts with just 1 mask and 2 passages, with great simplification and time savings. The varus-valgus angle is set by using the proper Varus-valgus Block (1), while the rotation is set by introducing a nail in the 0° or 3° (2) hole. Special inclinations may be obtained by rotating the mask manually, keeping the condilar prongs well attached to the condyles. Let the Varus-valgus Block slide through the guides, until both prongs touch the condyles. Fix the Orientation Mask (3) with nails. Once the mask is properly placed and fixed, let the Size Mask slide into the Orientation Mask (3), until the stylus touches the anterior part of the femur (4).

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

Verify that the size has not changed after the regulation of the stylus position. Introduce two headed nails into the hole corresponding to the chosen size. The device will automatically fit the cutting slots to the size.
 NOTE: in order to let the device find the best position, the nails must NOT be hammered, but pushed by hand, vibrating it a bit in the vertical direction to facilitate the sliding of the mechanisms. As the system has slid in the correct position, and the nail has reached the femoral bone, use a hammer to fix the nail to the bone. If the smaller size is chosen, the graduate scale will sign +2 or +4; these numbers will be used in the distal cut (step 14).

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Make the anterior and posterior cuts through the proper slots.

The stylus, when positioned at the correct size, signs the proximal border of the anterior cut.

Insert two headless nails in the upper holes, then remove all devices, leaving the nails in place. To remove the devices, remove all the Headed Nails, then lift up and remove the Size Mask, then extract the Orientation Mask and the Intramedullar Rod.

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The Distal Femoral Cutting Block presents 3 couples of holes, marked as 0, +2, +4.

It is possible to insert it in the proper couple of holes, previously determined (the number reported in the lower graduated scale (1) in figure A must coincide with the number of the hole).

NOTE: In any case, the final evaluation of the cut depth must be made during surgery by the surgeon, following the anatomy of the patient.
Cut the distal condyles by using the slot, then remove the Distal Cut Block, leaving the nails in place.

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At this stage, it is very convenient to perform the balance checks to verify the actual matching of the femoral and tibial parts. It is important to do it before chamfering, to obtain a more precise measurement and a much easier recutting, if any.




A very easy recutting is enabled by the presence of the nails in the original position.

Screwing the Male Control Rod on the Female Control Rod, check the Anatomic axis alignment.

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The Spacer Augmentations are marked with the thickness size (thickness reference number) of the correspondent Insert.
Read the thickness reference number on the used augmentation (+2, +4, +7). If the collateral ligament tension is satisfactory without using any augmentation, then use the +0 insert.

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
Make the finishing cut leaning the blade on the oblique plane, as showed in the figure to the left. In the case of a CR implant, make the holes for the pin and remove the mask.
In the case of the PS implant, make the box cuts leaning the blade on the oblique plane on the two lateral planes and the posterior plane, as showed in the figure to the right. Remove the mask, then refine the cuts with an osteotome or an axil blade, if necessary.
Don't make the holes for the pins, absent in the PS implant.

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Put a Trial Tibial Tray in place with the modular handle plugged. Change the tray size until you find the right one for the tibial component. Spike it with at least four headed nails. NOTE: The best sequence of insertion of the pins is: vertical posterior short nails; vertical anterior short nails; obliques short nails. Use the Checking Rod to check the actual alignment of the tibial cut. Put the Trial-Insert Adapter in its slot on the Trial Tibial Tray. Following the indications about the thickness and taking into account that, with the Mobile Hollow models, the femoral component must match the size of the insert exactly, while the tibial plate can be ± 1 size, choose the appropriate trial insert and place it on the Trial-Insert Adapter. Place the Trial Femoral Component and make a flexion and extension trial reduction to check the fitness of the whole implant.

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Remove all the trial components, except for the Trial Tibial Tray. Place the MBH Tibial Drill Guide (1) in its slot on the Trial Tibial Tray. Set the drill height (2) at the chosen size. Place the Arrow Broach in its slot on the Trial Tibial Tray. Set the Tibial Plate size on the upper graduate scale (2); hammer, then check the complete penetration of the broach with the lower graduate scale (3).

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After the definitive Tibial Plate is mounted, a further trial reduction may be performed using the Trial Insert again and the Tapered Adapter.

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Mount the final Femoral Component and the Insert, then reduce the articulation.

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If the implant of the Patellar Component is foreseen, the first step is to measure the thickness of the patella. Subtract 8-9-10 mm (the thickness of the three sizes of the patellar component) from the measured value, then set up the resulting value on the patellar pliers stylus for patella resection. By means of the Patellar Pliers, seize the patella, leaning it on the stylus. Pay attention to the inclination of the patellar resection. Cut the patella posterior part.

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Among the available patellar components, choose the one that better fits the resection. Holding the component by the side, make three holes with the Patellar Pin Drill.

Trial patellar components are available into the instrument set. During the cementation process, the Patellar Cementation Pliers can be used to keep the patellar component pressed onto the bone.

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