



## IMF Quick-Fix Set

The IMF Quick-Fix Set guarantees a quick, safe and stitch-free intermaxillary fixation in traumatology and orthognatic surgery. Operations in oral and maxillofacial surgery very often require a safe intermaxillary fixation in order to achieve a good and reliable occlusion. Up to now neither infections nor irritations of the mucosa have occurred during clinical use. Biocompatibility of the IMF screws and their atraumatic head allow a short as well as long-term IMF without problems.

Advantages:

- Protects surgeons from infection
- Time-saving
- Safe and reliable IMF even in the case of defective natural dentures or edentulous (toothless) patients
- Prevention against caries and injury to the gum, mostly in the case of long-term IMF
- Fast and easy removal
- Permanent availability, always ready to use





**IMF Screws** (Packaging: 4/ea.)



**MCD Screw, self-drilling**



**CROSS-FIT Screw, self-drilling**



Art. No.	Diameter x Length
25-20508	2.0 x 8 mm
25-20511	2.0 x 11 mm
25-20514	2.0 x 14 mm

Art. No.	Diameter x Length
25-20308	2.0 x 8 mm
25-20311	2.0 x 11 mm
25-20314	2.0 x 14 mm

**Instruments**



10-65115	Drill 1.5 x 50 mm, WL 22 mm, Stryker
----------	--------------------------------------



11-12026	MCD blade, self-retaining
----------	---------------------------



11-12025	CROSS-FIT blade, self-retaining
----------	---------------------------------



11-12000	Screwdriver handle, 10 cm
----------	---------------------------



25-09051	Stainless steel ligature wire, 0.5 mm, 10 m
----------	---



Container



73-31970 IMF tray incl. lid, length clips and case clip

IMF Quick-Fix Set (Recommendation) Ordering information

Art. No.	Quantity	Description <span style="float: right;">+</span>
<b>25-20305</b>	1	<b>IMF-Quick-Fix Set MCD</b> complete consisting of:
		<b>Sterilization and Storage</b>
73-31970	1	IMF tray incl. lid and length clips
		<b>Instruments</b>
11-12000	1	Screwdriver handle, 10 cm
11-12026	1	MCD blade
		<b>Drills and Accessories</b>
10-65115	1	Drill 1.5 x 50 mm, WL 22 mm, Stryker
25-09051	1	Stainless steel ligature wire 0.5 mm, 10 m
		<b>MCD Screws IMF</b>
25-20508	1	MCD Screw 2.0 x 8 mm, 4/ea.
25-20511	1	MCD Screw 2.0 x 11 mm, 4/ea.
25-20514	1	MCD Screw 2.0 x 14 mm, 4/ea.

Art. No.	Quantity	Description <span style="float: right;">+</span>
<b>25-20302</b>	1	<b>IMF-Quick-Fix Set CROSS-FIT</b> complete consisting of:
		<b>Sterilization and Storage</b>
73-31970	1	IMF tray incl. lid and length clips
		<b>Instruments</b>
11-12000	1	Screwdriver handle, 10 cm
11-12025	1	CROSS-FIT blade
		<b>Drills and Accessories</b>
10-65115	1	Drill 1.5 x 50 mm, WL 22 mm, Stryker
25-09051	1	Stainless steel ligature wire 0.5 mm, 10 m
		<b>CROSS-FIT Screws IMF</b>
25-20308	1	CF Screw 2.0 x 8 mm, 4/ea.
25-20311	1	CF Screw 2.0 x 11 mm, 4/ea.
25-20314	1	CF Screw 2.0 x 14 mm, 4/ea.

High Quality  
Implants  
Made by MONDEAL

We reserve the right to make catalog and design changes arising from further developments and modifications. Pictures, product descriptions and texts are the property of MONDEAL Medical Systems GmbH. Further use and reproduction by Third Parties requires written consent from MONDEAL Medical Systems GmbH. All rights reserved! Note: Responsibility for the selection of the patient, for adequate training as well as the decision as to the choice and use of the implants, a post-operative removal or retention of implants is the sole decision of the surgeon. All brands and logos referred to in this brochure, registered by Third Parties, are subject without restriction to the legal provisions currently in force and to the rights of ownership of the respective registered owner.

© 2011 - MONDEAL Medical Systems GmbH

