



# Logic JR.

PEDIATRIC DISTRACTION SYSTEM




### THE LOGICAL CURVILINEAR SOLUTION FOR PEDIATRIC APPLICATION

- THE PERFECT SOLUTION FOR NEONATES OR INFANTS UNDER 4 YEARS OF AGE TO CORRECT CONGENITAL DEFICIENCIES OR POST TRAUMATIC DEFECTS OF THE MANDIBLE
- Low profile intra-oral device for use in cases where distraction osteogenesis is indicated
- Utilizes a patented flexible wire that provides access through soft tissue for ease of activation
- Choice of two devices with fixed curvilinear trajectories that MIMIC THE NATURAL GROWTH CURVE OF THE MANDIBLE

PRODUCTS THAT IMPROVE PEOPLE'S LIVES

## LOGIC JR. PEDIATRIC DISTRACTION SYSTEM

Refer to Product Insert for specific indications and Instructions For Use.



#### **Implants**

Part Number	DESCRIPTION
216-1000	LOGIC™ Jr. Distractor, Left 52mm
216-1001	LOGIC™ Jr. Distractor, Right 52mm
216-1002	LOGIC™ Jr. Distractor, Straight (Left/Right)

### Instruments and Accessories

PART NUMBER	DESCRIPTION
216-1200	Logic™ Jr. Organizer Block
216-1201	0.9mm Pilot Drill, Logic™ Jr., OsteoPower
216-1202	0.9mm Pilot Drill, Logic™ Jr., J-Latch
220-0081	1.2mm Screwdriver Shaft, Long, Manual
216-0102	Distraction Tool
216-0103	Activation Wire Removal Tool
216-0305	Silicone Tubing Sleeve
216-0310	Logic/Logic Jr Template Kit
220-0019	TAPERLOCK™ Screwdriver Body
220-0022	1.2mm Screwdriver Shaft, Manual
220-0071	1.2mm Screw VISIDISK™ (YELLOW)
220-0085	1.2mm AUTO-DRIVE™ Screw VISIDISK™ (YELLOW)

#### Screws

_	
PART NUMBER	DESCRIPTION
211-1203	1.2mm x 3mm AUTO-DRIVE™ Screw
211-1204	1.2mm x 4mm AUTO-DRIVE™ Screw
211-1205	1.2mm x 5mm AUTO-DRIVE™ Screw
204-1503	1.5mm x 3mm Safety Screw
204-1506	1.5mm x 6mm Safety Screw

#### Activation Wires

Part Number	DESCRIPTION
216-0302	Activation Wire, 62mm
216-0306	Activation Wire, 72mm
216-0307	Activation Wire, 82mm



3885 Arapaho Road

Addison, Texas 75001

Customer Service: 800.456.7779

Main: 972.677.4600

Fax: 800.390.2620

www.osteomedcorp.com

P/N 030-1439 Rev. B

